


Anti-Samm50 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Y512

Item Number: CM0004966



Introduction

High-quality rabbit monoclonal antibody targeting human SAMM50, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	SAMM50 Rabbit mAb
Remarks/Alias	OMP85; SAM50; TOB55; TRG-3; CGI-51; YNL026W; SAMM50
Species	Human
GeneID	25813
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human SAMM50 (Q9Y512).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y512
Protein Weight	52kDa
Shipping	Ice bag