

Anti-Sap130 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H0E3

Item Number: CM0017889

Introduction

Rabbit polyclonal antibody against human SAP130, validated for WB and ELISA with cross-reactivity to mouse. Ideal for studying transcriptional repression and mSin3A complex function.

[Learn More](#)

Parameter	Value
Product Name	SAP130 Rabbit pAb
Remarks/Alias	SAP130
Species	Human
Gene ID (Human)	79595
Gene ID	79595
Immunogen	Recombinant protein/Recombinant fusion protein containing a sequence corresponding to amino acids 819-1048 of human SAP130 (NP_078821.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9H0E3
Protein Weight	110kDa
Shipping	Ice bag