


Anti-Sar1B Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y6B6

Item Number: CM0030917



Introduction

Rabbit polyclonal antibody against human SAR1B (SWISS: Q9Y6B6). Validated for WB and ELISA. Targets the small COPII coat GTPase involved in ER-to-Golgi transport and lipoprotein trafficking. Ideal for diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	SAR1B Rabbit pAb
Remarks/Alias	ANDD; CMRD; GTBPB; SARA2; SAR1B
Species	Human
GeneID (Human)	51128
GeneID	51128
Immunogen	Recombinant fusion protein corresponding to amino acids 1-198 of human SAR1B (NP_057187.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y6B6
Protein Weight	22kDa
Shipping	Ice bag