

Anti-Sarcospan (Sspn) Rabbit Polyclonal Antibody For Wb - Q14714

Item Number: CM0013698

Introduction

Rabbit polyclonal antibody against human Sarcospan (SSPN), validated for WB, ELISA. Cross-reacts with mouse, rat. Suitable for dystrophin-glycoprotein complex and muscle biology research.

[Learn More](#)

Parameter	Value
Product Name	SSPN Rabbit pAb
Remarks/Alias	KRAG; NSPN; SPN1; SPN2; DAGA5; SSPN
Species	Human
GeneID	8082
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-60 of human SSPN (NP_005077.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14714
Protein Weight	27kDa
Shipping	Ice bag