

Anti-Sardh Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9UI12

Item Number: CM0028925

Introduction

Rabbit polyclonal antibody against human SARDH (sarcosine dehydrogenase). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for mitochondrial choline metabolism research.

[Learn More](#)

Parameter	Value
Product Name	SARDH Rabbit pAb
Remarks/Alias	SAR; SDH; SARD; BPR-2; DMGDHL1; SARDH
Species	Human
Gene ID (Human)	1757
Gene ID	1757
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-560 of human SARDH (NP_009032.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UL12
Protein Weight	101kDa
Shipping	Ice bag