

Anti-Sarm1 Monoclonal Antibody For Wb, Elisa - Q6Pds3

Item Number: CM0007083

Introduction

Mouse monoclonal antibody targeting Sarm1, validated for WB and ELISA. Detects endogenous Sarm1 in human, mouse, and rat samples. Suitable for axonal degeneration and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	Sarm1 Mouse mAb
Remarks/Alias	Sarm; MyD885; A830091I15Rik; Sarm1
Species	Mouse
Gene ID (Human)	23098
Gene ID	237868
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse Sarm1 (NP_766383.2).
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6PDS3
Protein Weight	79kDa
Shipping	Ice bag