


Anti-Sarm1 Polyclonal Antibody For Wb, Ip, Elisa - Q6Szw1

Item Number: CM0022419



Introduction

High-quality rabbit polyclonal anti-SARM1 antibody for Western blot, IP, and ELISA. Detects human SARM1, a key NAD⁺ hydrolase in programmed axon death and innate immune signaling.

[Learn More](#)

Parameter	Value
Product Name	SARM1 Rabbit pAb
Remarks/Alias	SARM; HsTIR; SAMD2; hSARM1; MyD88-5; SARM1
Species	Human
Gene ID (Human)	23098
Gene ID	23098
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 625-724 of human SARM1 (NP_055892.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q6SZW1
Protein Weight	79kDa
Shipping	Ice bag