


Anti-Sarnp Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P82979

Item Number: CM0015228



Introduction

Rabbit polyclonal antibody against human SARNP (CIP29/Hcc-1), validated for WB, IHC-P, IF/ICC, and ELISA. SARNP is a TREX complex component involved in mRNA export and cell growth regulation. Suitable for cancer and transcriptional research.

[Learn More](#)

Parameter	Value
Product Name	SARNP Rabbit pAb
Remarks/Alias	HCC1; THO1; CIP29; HSPC316; SARNP
Species	Human
Gene ID (Human)	84324
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human SARNP (NP_149073.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P82979
Protein Weight	24kDa
Shipping	Ice bag