

Anti-Sars-Cov-2 3CLpro Monoclonal Antibody For Wb, If/Icc, Elisa - P0Dtc2

Item Number: CM0004555



Introduction

SARS-CoV-2 3CLpro rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA. Recognizes SARS-CoV-2 3CLpro protein. Recombinant immunogen corresponding to amino acids 1-306 of 3CLpro. Ideal for COVID-19 research.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 3CLpro Rabbit mAb
Gene ID	43740578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-306 of coronavirus SARS-CoV-2 3CLpro (YP_009725295.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTC2
Protein Weight	141kDa
Shipping	Ice bag