

Anti-Sars-Cov-2 3Clpro Polyclonal Antibody For Wb, If/Icc, Elisa - P0Dtc2

Item Number: CM0024409



Introduction

High-quality rabbit polyclonal antibody against SARS-CoV-2 3CLpro, the main viral protease. Ideal for WB, IF/ICC, and ELISA applications in COVID-19 research, therapeutic monitoring, and diagnostic kit development.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 3CLpro Rabbit pAb
GeneID	43740578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-306 of coronavirus 3CLpro (YP_009725295.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTC2
Protein Weight	141kDa
Shipping	Ice bag