


Anti-Sars-Cov-2 N Protein Monoclonal Antibody For Wb, Elisa - P59595

Item Number: CM0003192



Introduction

High-quality rabbit monoclonal antibody against SARS-CoV-2 N Protein (Nucleoprotein), validated for WB and ELISA. Ideal for COVID-19 research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 N Protein Rabbit mAb, BSA and glycerol free
Remarks/Alias	Nucleoprotein; NP; SARS-CoV-2 N Protein
Species	coronavirus
GenelD (human)	2019-nCoV
GenelD	43740575
Immunogen	Recombinant protein; Recombinant fusion protein of SARS-CoV-2 Nucleoprotein Full length (P59595)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P59595
Protein Weight	46kDa
Shipping	Dry ice