


Anti-Sars-Cov-2 N Protein Monoclonal Antibody For Wb, If/Icc, Elisa - P59595

Item Number: CM0002945



Introduction

SARS-CoV-2 N Protein Mouse mAb validated for WB, IF/ICC, ELISA. Recognizes SARS-CoV-2 nucleoprotein with cross-reactivity to SARS-CoV. Suitable for virology, diagnostics, and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 N Protein Mouse mAb
Remarks/Alias	Nucleoprotein; NP; SARS-CoV-2 N Protein
Species	coronavirus
Gene ID (Human)	SARS-CoV-2
Gene ID	43740575
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-420 of coronavirus Nucleoprotein (P0DTC9)
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P59595
Protein Weight	46kDa
Shipping	Ice bag