


## Anti-Sars-Cov-2 Nsp12 Rabbit Polyclonal Antibody For Wb, Elisa - P0Dtd1

Item Number: CM0020123



### Introduction

High-quality rabbit polyclonal antibody against SARS-CoV-2 NSP12 (RNA polymerase) for Western blot (WB) and ELISA. Recognizes SARS-CoV-2; useful for studying viral replication and diagnostic development.

[Learn More](#)

| Parameter        | Value                                                                                                          |
|------------------|----------------------------------------------------------------------------------------------------------------|
| Product Name     | SARS-CoV-2 NSP12 Rabbit pAb                                                                                    |
| GeneID           | 43740578                                                                                                       |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 200-300 of coronavirus NSP12 (YP_009725307.1) |
| Source           | Rabbit                                                                                                         |
| Category         | Polyclonal Antibodies                                                                                          |
| Application      | WB, ELISA                                                                                                      |
| Cross Reactivity | SARS-CoV-2                                                                                                     |
| SWISS            | P0Dtd1                                                                                                         |
| Protein Weight   | 794kDa                                                                                                         |
| Shipping         | Ice bag                                                                                                        |