


Anti-Sars-Cov-2 Nsp13 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P0Dtd1

Item Number: CM0022063



Introduction

Rabbit polyclonal antibody against SARS-CoV-2 NSP13 helicase, validated for WB, IF/ICC, ELISA. Cross-reacts with SARS-CoV-2. Ideal for viral replication and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 NSP13 Rabbit pAb
Species	SARS-CoV-2
GeneID	43740578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-601 of coronavirus NSP13 (YP_009725308.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0Dtd1
Protein Weight	794kDa
Shipping	Ice bag