


Anti-Sars-Cov-2 Nsp3 Polyclonal Antibody For Wb, Elisa - P0Dtd1

Item Number: CM0020741



Introduction

SARS-CoV-2 NSP3 rabbit polyclonal antibody validated for WB and ELISA. Targets the NSP3 region of the SARS-CoV-2 replicase polyprotein 1ab (UniProt P0DTD1). Ideal for viral replication and host interaction studies.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 NSP3 Rabbit pAb
GeneID	43740578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of coronavirus NSP3 (YP_009725299.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTD1
Protein Weight	794kDa
Shipping	Ice bag