


Anti-Sars-Cov-2 Nsp4 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P0Dtd1

Item Number: CM0022837



Introduction

Anti-SARS-CoV-2 NSP4 rabbit polyclonal antibody suitable for WB, IP, ELISA. Targets a key viral replicase protein involved in double-membrane vesicle formation, aiding SARS-CoV-2 research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 NSP4 Rabbit pAb
Gene ID	43740578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-500 of coronavirus NSP4 (YP_009725300.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTD1
Protein Weight	794 kDa
Shipping	Ice bag