

Anti-Sars-Cov-2 Nsp7 Rabbit Polyclonal Antibody For Wb, Elisa - P0Dtc2

Item Number: CM0021254



Introduction

Rabbit polyclonal antibody against SARS-CoV-2 NSP7 protein. Validated for Western blot (WB) and ELISA, with human cross reactivity. Suitable for COVID-19 research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 NSP7 Rabbit pAb
Gene ID	43740578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-83 of coronavirus NSP7 (YP_009742614.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P0DTC2
Protein Weight	141kDa
Shipping	Ice bag