


Anti-Sars-Cov-2 Nsp9 Rabbit Polyclonal Antibody For Wb, Elisa - P0Dtd1

Item Number: CM0020091



Introduction

Rabbit polyclonal antibody raised against recombinant SARS-CoV-2 NSP9 (replicase polyprotein 1ab, amino acids 1-113). Suitable for Western blot and ELISA applications. Recognizes the NSP9 region of the viral replicase polyprotein (Swiss-Prot P0DTD1).

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 NSP9 Rabbit pAb
GeneID	43740578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of coronavirus NSP9 (YP_009742616.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTD1
Protein Weight	794kDa
Shipping	Ice bag