

Anti-Sars-Cov-2 Nucleoprotein Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P0Dtc9

Item Number: CM0003772



Introduction

Rabbit monoclonal antibody against SARS-CoV-2 nucleoprotein (N protein), validated for WB, IF/ICC, and ELISA. Suitable for COVID-19 research, diagnostic assay development, and vaccine studies. Target accession: P0DTC9.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 Nucleoprotein Rabbit mAb
Remarks/Alias	N; Nucleoprotein; NC; SARS-CoV-2 Nucleoprotein
Species	Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2)
GeneID	43740575
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTC9
Protein Weight	46kDa
Shipping	Ice bag