


Anti-Sars-Cov-2 Orf7A Polyclonal Antibody For Wb And Elisa - P0Dtc7

Item Number: CM0021405



Introduction

SARS-CoV-2 ORF7a rabbit polyclonal antibody validated for WB and ELISA. Recognizes ORF7a (P0DTC7), a viral immunoevasin that inhibits host immune responses by antagonizing tetherin and downregulating MHC-I. Ideal for COVID-19 research.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 ORF7a Rabbit pAb
Species	SARS-CoV-2
Gene ID	43740573
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-98 of coronavirus ORF7a (YP_009724395.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTC7
Protein Weight	14kDa
Shipping	Ice bag