

Anti-Sars-Cov-2 Spike Ecd Polyclonal Antibody For Wb, Elisa - P0Dtc2

Item Number: CM0021544



Introduction

SARS-CoV-2 Spike ECD rabbit polyclonal antibody validated for Western blot and ELISA. Recognizes the extracellular domain of spike glycoprotein, essential for viral entry research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 Spike ECD Rabbit pAb
Remarks/Alias	spike glycoprotein; SARS-CoV-2 Spike ECD
GeneID	43740568
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 11-1208 of coronavirus Spike ECD (YP_009724390.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTC2
Protein Weight	141kDa
Shipping	Ice bag