

Anti-Sars-Cov-2 Spike Rabbit Polyclonal Antibody For Wb - P0Dtc2

Item Number: CM0020686



Introduction

Rabbit polyclonal antibody targeting SARS-CoV-2 Spike glycoprotein (S). Validated for WB, DB, IF/ICC, IP, and ELISA. Specific for SARS-CoV-2 Spike. Suitable for viral entry, fusion, and immunoevasion research.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 Spike Rabbit pAb
Remarks/Alias	spike glycoprotein; SARS-CoV-2 Spike
GeneID	43740568
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 600-700 of coronavirus Spike (YP_009724390.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, DB, IF/ICC, IP, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTC2
Protein Weight	141kDa
Shipping	Ice bag