

Anti-Sars-Cov-2 Spike Rbd (Omicron) Polyclonal Antibody For Wb, Elisa - P0Dtc2

Item Number: CM0020630



Introduction

Rabbit polyclonal antibody targeting the receptor-binding domain of SARS-CoV-2 Spike protein (Omicron variant). Suitable for Western blot and ELISA. Recognizes SARS-CoV-2. Ideal for COVID-19 research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	SARS-COV-2 Spike RBD(Omicron) Protein pAb
Remarks/Alias	spike glycoprotein
Gene ID	43740568
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 319-541 of coronavirus SARS-COV-2 Spike RBD(Omicron) Protein (P0DTC2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTC2
Protein Weight	141kDa
Shipping	Ice bag