

Anti-Sars-Cov-2 Spike Rbd Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P0Dtc2

Item Number: CM0023226



Introduction

Rabbit polyclonal antibody recognizing the receptor-binding domain (RBD) of SARS-CoV-2 spike glycoprotein (P0DTC2). Validated for WB, IF/ICC, IP, and ELISA; reacts with SARS-CoV-2. Ideal for viral entry and neutralization studies.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 Spike RBD Rabbit pAb
Remarks/Alias	spike glycoprotein; SARS-CoV-2 Spike RBD
Gene ID	43740568
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 332-550 of coronavirus Spike RBD (YP_009724390.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTC2
Protein Weight	141kDa
Shipping	Ice bag