

Anti-Sars-Cov-2 Spike S1 Rabbit Polyclonal Antibody For Wb, Elisa - P59594

Item Number: CM0021770



Introduction

Rabbit polyclonal antibody against SARS-CoV-2 Spike S1 (E2), suitable for WB and ELISA. Cross-reacts with SARS-CoV. Target protein is 139 kDa; UniProt Accession: P59594.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 Spike S1 Rabbit pAb
Remarks/Alias	E2; SARS-CoV-2 Spike S1
Gene ID	1489668
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of coronavirus Spike S1 (NP_828851.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV
SWISS	P59594
Protein Weight	139kDa
Shipping	Ice bag