

Anti-Sars-Cov-2 Spike S1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P0Dtc2

Item Number: CM0021196



Introduction

This rabbit polyclonal antibody targets the SARS-CoV-2 Spike S1 subunit (UniProt P0DTC2), validated for WB, IF/ICC, IP, and ELISA. Ideal for viral entry and host interaction studies.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 Spike S1 Rabbit pAb
Remarks/Alias	spike glycoprotein; SARS-CoV-2 Spike S1
Gene ID	43740568
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 11-682 of coronavirus Spike S1 (YP_009724390.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, SARS-CoV-2
SWISS	P0DTC2
Protein Weight	141kDa
Shipping	Ice bag