

Anti-Sars-Cov-2 Spike S2 Ecd Rabbit Polyclonal Antibody For Wb, Db, If/Icc, Ip, And Elisa - P0Dtc2

Item Number: CM0023544



Introduction

Rabbit polyclonal antibody targeting SARS-CoV-2 Spike S2 extracellular domain (ECD). Validated for WB, DB, IF/ICC, IP, and ELISA. Recognizes human and SARS-CoV-2. Suitable for viral entry and fusion research.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 Spike S2 ECD Rabbit pAb
Remarks/Alias	spike glycoprotein; SARS-CoV-2 Spike S2 ECD
Species	Rabbit
Gene ID	43740568
Immunogen	Synthetic peptide corresponding to amino acids 686-1214 of SARS-CoV-2 Spike S2 ECD (YP_009724390.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, DB, IF/ICC, IP, ELISA
Cross Reactivity	Human, SARS-CoV-2
SWISS	P0DTC2
Protein Weight	141kDa
Shipping	Ice bag