

Anti-Sars-Cov-2 Spike S2 Polyclonal Antibody For Wb, Ip, Elisa - P0Dtc2

Item Number: CM0024350



Introduction

Rabbit polyclonal antibody against SARS-CoV-2 Spike S2 subunit, suitable for Western blot, immunoprecipitation, and ELISA. Detects spike glycoprotein of SARS-CoV-2 (UniProt P0DTC2) with cross-reactivity to SARS-CoV-2.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 Spike S2 Rabbit pAb
Remarks/Alias	spike glycoprotein; SARS-CoV-2 Spike S2
Gene ID	43740568
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTC2
Protein Weight	141kDa
Shipping	Ice bag