

Anti-Sars-Cov Spike Rbd Monoclonal Antibody For Wb, Elisa - P59594

Item Number: CM0006970



Introduction

Recombinant rabbit monoclonal antibody targeting the receptor-binding domain (RBD) of SARS-CoV Spike glycoprotein (P59594). Validated for Western blot and ELISA applications. Suitable for viral entry and immunoevasion research.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV Spike RBD Rabbit mAb
Remarks/Alias	E2; SARS-CoV Spike RBD
Gene ID	1489668
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 350-450 of coronavirus SARS-CoV Spike RBD (NP_828851.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV
SWISS	P59594
Protein Weight	139kDa
Shipping	Ice bag