


Anti-Sars Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P49591

Item Number: CM0016249



Introduction

Rabbit polyclonal antibody against human SARS (seryl-tRNA synthetase), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for studying translation and angiogenesis regulation.

[Learn More](#)

Parameter	Specification
Product Name	SARS Rabbit pAb
Remarks/Alias	SARS; SERS; SERRS; NEDMAS
Species	Human
Gene ID (Human)	6301
Gene ID	6301
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SARS (NP_006504.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49591
Protein Weight	59kDa
Shipping	Ice bag