


Anti-Sars/Seryl-Trna Synthetase (Sars1) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49591

Item Number: CM0029703



Introduction

High-quality rabbit polyclonal antibody against SARS (seryl-tRNA synthetase, SARS1) for Western blot, IHC-P, IF/ICC, and ELISA. Shows cross-reactivity with human, mouse, rat. Essential for studying translational and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	SARS Rabbit pAb
Remarks/Alias	SARS; SERS; SERRS; NEDMAS
Species	Human
GeneID (Human)	6301
GeneID	6301
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SARS (NP_006504.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49591
Protein Weight	59kDa
Shipping	Ice bag