


Anti-Sart1 Polyclonal Antibody (Rabbit Pab) For Wb, Ihc-P, Ip, Elisa - O43290

Item Number: CM0029059



Introduction

Rabbit polyclonal antibody against human SART1, a U4/U6.U5 tri-snRNP splicing factor. Validated for WB, IHC-P, IP, and ELISA applications. Cross-reacts with mouse and rat. Suitable for diagnostic and research use.

[Learn More](#)

Parameter	Value
Product Name	SART1 Rabbit pAb
Remarks/Alias	HAF; Ara1; HOMS1; Snu66; SART1259; SNRNP110; SART1
Species	Human
GeneID	9092
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43290
Protein Weight	90kDa
Shipping	Ice bag