


Anti-Sart3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15020

Item Number: CM0014680



Introduction

Rabbit polyclonal antibody against SART3, a spliceosome recycling factor. Validated for WB, IHC-P, ELISA in human, mouse, rat. Ideal for studying mRNA splicing and gene regulation.

[Learn More](#)

Parameter	Value
Product Name	SART3 Rabbit pAb
Remarks/Alias	P100; p110; DSAP1; TIP110; p110(nrb); RP11-13G14; SART3
Species	Human
Gene ID (Human)	9733
Gene ID	9733
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 664-963 of human SART3 (NP_055521.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15020
Protein Weight	110kDa
Shipping	Ice bag