


Anti-Sat1 Rabbit Polyclonal Antibody For Wb, Elisa - P21673

Item Number: CM0028223



Introduction

High-quality rabbit polyclonal antibody against human SAT1, validated for Western blot and ELISA. Detects SAT1 also known as SSAT, involved in polyamine acetylation. Cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	SAT1 Rabbit pAb
Remarks/Alias	SAT, DC21, KFSD, SSAT, KFSDX, SSAT-1, SAT1
Species	Human
GeneID (Human)	6303
GeneID	6303
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-171 of human SAT1 (NP_002961.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P21673
Protein Weight	20kDa
Shipping	Ice bag