


Anti-Satb1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q01826

Item Number: CM0006660



Introduction

High-specificity rabbit monoclonal antibody against SATB1, suitable for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Cross-reacts with human, mouse, and rat SATB1. For chromatin and transcription research.

[Learn More](#)

Parameter	Value
Product Name	SATB1 Rabbit mAb
Remarks/Alias	DEFDA; KTZSL; SATB1
Species	Human
GeneID (Human)	6304
GeneID	6304
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 572-763 of human SATB1 (NP_002962.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01826
Protein Weight	86kDa
Shipping	Ice bag