


Anti-Satb1 Rabbit Polyclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - Q01826

Item Number: CM0024767



Introduction

SATB1 rabbit polyclonal antibody validated for WB, IP, IF-P, IHC-P, ELISA. Reacts with human, mouse, and rat. Suitable for chromatin remodeling, T-cell development, and cancer metastasis studies.

[Learn More](#)

Parameter	Value
Product Name	SATB1 Rabbit pAb
Remarks/Alias	DEFDA; KTZSL; SATB1
Species	Human
Gene ID (Human)	6304
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 572-763 of human SATB1 (NP_002962.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01826
Protein Weight	86 kDa
Shipping	Ice bag