


Anti-Satb2 Monoclonal Antibody For Wb, If-P, Ihc-P, Ip, Elisa - Q9Upw6

Item Number: CM0003591



Introduction

Rabbit monoclonal antibody against human SATB2, validated for WB, IF, IHC, IP, and ELISA. Cross-reacts with mouse. Ideal for chromatin and neurodevelopment studies.

[Learn More](#)

Parameter	Value
Product Name	SATB2 Rabbit mAb
Remarks/Alias	GLSS, DEL2Q32Q33, C2DELq32q33, SATB2
Species	Human
GeneID (human)	23314
GeneID	23314
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 634-733 of human SATB2 (NP_001165980.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UPW6
Protein Weight	83kDa
Shipping	Ice bag