


Anti-Satb2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9Upw6

Item Number: CM0003124



Introduction

High-performance anti-SATB2 rabbit monoclonal antibody validated for WB, IHC-P, IP, and ELISA. Reacts with human, mouse, and rat SATB2, a key transcription factor in neuronal and skeletal development.

[Learn More](#)

Parameter	Value
Product Name	SATB2 Rabbit mAb
Remarks/Alias	GLSS; DEL2Q32Q33; C2DELq32q33; SATB2
Species	Human
GeneID (Human)	23314
GeneID	23314
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human SATB2 (Q9UPW6).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UPW6
Protein Weight	83 kDa
Shipping	Ice bag