


Anti-Satb2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Upw6

Item Number: CM0022425



Introduction

High-quality rabbit polyclonal antibody against human SATB2 (SWISS: Q9UPW6). Validated for Western blot, immunohistochemistry (paraffin), and ELISA. Suitable for research on transcription regulation, neurodevelopment, and osteoblast differentiation.

[Learn More](#)

Parameter	Value
Product Name	SATB2 Rabbit pAb
Remarks/Alias	GLSS; DEL2Q32Q33; C2DELq32q33; SATB2
Species	Human
Gene ID (Human)	23314
Gene ID	23314
Immunogen	Synthetic peptide corresponding to amino acids 634-733 of human SATB2 (NP_001165980.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9UPW6
Protein Weight	83kDa
Shipping	Ice bag