

Anti-Sav1 Polyclonal Antibody For Wb, Elisa - Q9H4B6

Item Number: CM0029977



Introduction

Rabbit polyclonal antibody against human SAV1 (hWW45), validated for WB and ELISA. Cross-reacts with mouse. Ideal for Hippo signaling, organ size control, and tumor suppression research.

[Learn More](#)

Parameter	Value
Product Name	SAV1 Rabbit pAb
Remarks/Alias	SAV; WW45; WWP4; SAV1
Species	Human
GeneID (human)	60485
GeneID	60485
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 199-383 of human SAV1 (NP_068590.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H4B6
Protein Weight	45kDa
Shipping	Ice bag