

Anti-Sbf1 Rabbit Polyclonal Antibody For Wb, Elisa - 095248

Item Number: CM0013758

Introduction

Rabbit polyclonal antibody targeting human SBF1 (MTMR5), validated for WB and ELISA applications, with cross-reactivity to human and rat. Ideal for investigating autophagy, neuronal myelination, and Charcot-Marie-Tooth disease pathways.

[Learn More](#)

Parameter	Value
Product Name	SBF1 Rabbit pAb
Remarks/Alias	MTMR5; CMT4B3; DENND7A; SBF1
Species	Human
GeneID (human)	6305
GeneID	6305
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1630-1770 of human SBF1 (NP_002963.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O95248
Protein Weight	208kDa
Shipping	Ice bag