

Anti-Scaf4 Rabbit Polyclonal Antibody For Ihc-P And Elisa - O95104

Item Number: CM0023108



Introduction

High-quality rabbit anti-SCAF4 polyclonal antibody validated for IHC-P and ELISA. Recognizes human, mouse, and rat SCAF4. Ideal for studying transcriptional regulation and premature mRNA termination mechanisms.

[Learn More](#)

Parameter	Value
Product Name	SCAF4 Rabbit pAb
Remarks/Alias	SRA4; SFRS15; SCAF4
Species	Human
Gene ID (Human)	57466
Gene ID	57466
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-248 of human SCAF4 (NP_065757.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95104
Protein Weight	126kDa
Shipping	Ice bag