


Anti-Scaf8 (Rbm16) Rabbit Polyclonal Antibody For Wb, Elisa - Q9Upn6

Item Number: CM0021441



Introduction

Rabbit polyclonal antibody targeting human SCAF8 (RBM16), suitable for Western blotting and ELISA. Cross-reacts with mouse and rat. Useful for transcription termination and RNA processing studies. Recombinant immunogen, 141 kDa.

[Learn More](#)

Parameter	Value
Product Name	SCAF8 Rabbit pAb
Remarks/Alias	RBM16; SCAF8
Species	Human
Gene ID (Human)	22828
Gene ID	22828
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 572-634 of human SCAF8 (NP_055707.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UPN6
Protein Weight	141kDa
Shipping	Ice bag