

Anti-Scamp1 Rabbit Polyclonal Antibody For Wb, Elisa - O15126

Item Number: CM0027770

Introduction

Rabbit polyclonal antibody against SCAMP1, validated for WB and ELISA. Recognizes human, mouse, rat SCAMP1. Ideal for studying post-Golgi recycling and membrane trafficking pathways.

[Learn More](#)

Parameter	Value
Product Name	SCAMP1 Rabbit pAb
Remarks/Alias	SCAMP; SCAMP37; SCAMP1
Species	Human
Gene ID (human)	9522
Gene ID	9522
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human SCAMP1 (NP_004857.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15126
Protein Weight	38kDa
Shipping	Ice bag