

Anti-Scamp2 Polyclonal Antibody For Wb, Elisa - 015127

Item Number: CM0028563

Introduction

Anti-SCAMP2 rabbit polyclonal antibody raised against recombinant human SCAMP2. Validated for Western blot and ELISA, detects human SCAMP2 (~37kDa) involved in post-Golgi recycling pathways.

[Learn More](#)

Product Name	SCAMP2 Rabbit pAb
Remarks/Alias	SCAMP2
Species	Human
Gene ID (Human)	10066
Gene ID	10066
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-152 of human SCAMP2 (NP_005688.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15127
Protein Weight	37kDa
Shipping	Ice bag