

Anti-Scand1 Rabbit Polyclonal Antibody For Wb, Elisa - P57086

Item Number: CM0029478

Introduction

High-quality anti-SCAND1 rabbit polyclonal antibody for Western blot and ELISA applications. Validated in human, mouse, and rat samples. Recombinant immunogen. For transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	SCAND1 Rabbit pAb
Remarks/Alias	RAZ1; SDP1; SCAND1
Species	Human
Gene ID (Human)	51282
Gene ID	51282
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-179 of human SCAND1 (NP_057642.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P57086
Protein Weight	19kDa
Shipping	Ice bag