

Anti-Sccpdh Rabbit Polyclonal Antibody For Wb, Elisa - Q8Nbx0

Item Number: CM0021662

Introduction

Rabbit polyclonal antibody targeting human SCCPDH (NET11, CGI-49), validated for WB and ELISA. Reacts with human and mouse; immunogen from recombinant protein. Ideal for studying oxidoreductase function.

[Learn More](#)

Parameter	Value
Product Name	SCCPDH Rabbit pAb
Remarks/Alias	NET11; CGI-49; SCCPDH
Species	Human
GeneID (Human)	51097
GeneID	51097
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human SCCPDH (NP_057086.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NBX0
Protein Weight	47kDa
Shipping	Ice bag