


## Anti-Scd Polyclonal Antibody For Wb, If/Icc, Elisa - 000767

Item Number: CM0024911



### Introduction

Rabbit polyclonal antibody against human Stearoyl-CoA desaturase (SCD). Validated for WB, IF/ICC, ELISA. Reacts with human and rat. Suitable for lipid metabolism and fatty acid biosynthesis research.

[Learn More](#)

| Parameter        | Value                                   |
|------------------|-----------------------------------------|
| Product Name     | [KD Validated] SCD Rabbit pAb           |
| Remarks/Alias    | SCD1; FADS5; SCDOS; hSCD1; MSTP008; SCD |
| Species          | Human                                   |
| GeneID (human)   | 6319                                    |
| GeneID           | 6319                                    |
| Immunogen        | synthetic peptide                       |
| Source           | Rabbit                                  |
| Category         | Polyclonal Antibodies                   |
| Application      | WB, IF/ICC, ELISA                       |
| Cross Reactivity | Human, Rat                              |
| SWISS            | 000767                                  |
| Protein Weight   | 42kDa                                   |
| Shipping         | Ice bag                                 |