


Anti-Scd1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - 000767

Item Number: CM0008006



Introduction

Rabbit monoclonal antibody against human SCD1 (stearoyl-CoA desaturase), validated for WB, IF/ICC, and ELISA. Targets fatty acid biosynthesis enzyme; ideal for metabolic and lipid research. Human cross-reactivity.

[Learn More](#)

Parameter	Value
Product Name	SCD1 Rabbit mAb
Remarks/Alias	SCD1; FADS5; SCDOS; hSCD1; MSTP008
Species	Human
Gene ID (Human)	6319
Gene ID	6319
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 260-359 of human SCD1 (NP_005054.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O00767
Protein Weight	42kDa
Shipping	Ice bag