


Anti-Scd1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00767

Item Number: CM0020711



Introduction

Rabbit polyclonal antibody targeting human SCD1, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Stearoyl-CoA desaturase involved in lipid biosynthesis and energy homeostasis.

[Learn More](#)

Parameter	Value
Product Name	SCD1 Rabbit pAb
Remarks/Alias	SCD1; FAD55; SCDOS; hSCD1; MSTP008; SCD
Species	Human
GeneID (human)	6319
GeneID	6319
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SCD (NP_005054.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00767
Protein Weight	42kDa
Shipping	Ice bag